



# Compatibilidad del MED

Q1 2026

v1.0.0

# El futuro de la implantología ya está aquí

## Descripción del producto y configuraciones disponibles

El Dispositivo Electromagnético Miniaturizado (MED) de Magdent tiene la misma forma y se coloca igual que un pilar de cicatrización convencional. En su interior, libera un campo electromagnético pulsado (PEMF) controlado durante 30 días.

Esta tecnología, científicamente validada, estimula, acelera y mejora la osteogénesis (formación ósea) hasta en un 60% y favorece la osteointegración (conexión hueso-implante) hasta en un 50%, contribuyendo a reducir el período de cicatrización hasta en dos tercios tras la colocación del implante y permitiendo la colocación de la corona en un tercio del tiempo habitual, según estudios publicados y revisados por pares.

Lo que Magdent modifica es el entorno biológico de cicatrización inmediatamente después de la colocación del implante. Al influir de forma activa en la inflamación temprana, la angiogénesis y la actividad de los osteoblastos, la respuesta tisular se vuelve más consistente, incluso en hueso comprometido.

El MED se inserta en el implante como un pilar de cicatrización convencional, sin modificar el protocolo quirúrgico, y es compatible con el 80% de las marcas y modelos de implantes.



# Lista de compatibilidad

Implant brand	System/model	Connection	Platform Type	Implant Diameter $\emptyset$	MED Nr.
AB Dental	I Types systems	Internal Hex Connection	Narrow Platform	3.75	MED257
AB Dental	I Types systems	Internal Hex Connection	Wide Platform	4.5	MED387
ACE	Infinity Octagon		Regular Platform 4.8	3.3	MED363
ACE	Infinity Octagon		Regular Platform 4.8	4.1	MED363
ACE	Infinity Octagon		Regular Platform 4.8	4.8	MED363
ACE	Infinity	Internal Hex Connection	Standard Platform	3.7 / 4.1	MED257
ACE	Infinity	Internal Hex Connection	Wide Platform	4.7 / 5.1	MED387
ACE	Tri-Cam (like Nobel)	Conical Connection	Narrow Platform	3.5	MED271
ACE	Tri-Cam (like Nobel)	Conical Connection		4.3	MED295
ACE	Tri-Cam (like Nobel)	Conical Connection		5	MED301
Adin	Touareg	Internal Hex Connection		3.5 / 3.75 / 4.2 / 5	MED257
Adin	Swell	Internal Hex Connection		3.5 / 3.75 / 4.2 / 5	MED257
Adin	Touareg close fit	Conical Connection	Wide Platform	4.3 / 5	MED356
AlfaGate	AlfaGate Bioactive	Internal Hex Connection	Standard Platform	3.3 / 3.75 / 4.2 / 4.7 / 5.2	MED257
AlfaGate	AlfaGate Bioactive	Internal Hex Connection	Standard Platform	6	MED257
Alfa-Gate	Conical		Regular Platform	4.3/5	MED356
Alpha Bio	SPI,DFI,ATID	Internal Hex Connection		3.3 / 3.75 / 4.2 / 5	MED257
Alpha Bio	ICE	Internal Hex Connection	Narrow Platform	3.7	MED257
Alpha Bio	ICE	Internal Hex Connection		3.75 / 4.2 / 4.6 / 5.3	MED257
Alpha Dent	Bio	Internal Hex Connection		3.75 / 4.2 / 5	MED257
Alpha Dent	Classic	Internal Hex Connection		3.75 / 4.2 / 5	MED257
Alpha Dent	Active	Internal Hex Connection		3.75 / 4.2 / 5	MED257
Alpha Dent	Conus	Internal Hex Connection		3.75 / 4.2 / 5	MED257
Ancladen	Anclalock		3.5	3.75/4.25/5	MED257
Anthogyr	Ossfit		4.8	3.5/4.2	MED363
ARDS	CIT	Internal Hex Connection	Standard Platform	3.75 / 4.2 / 5	MED257
ARDS	Smart	Internal Hex Connection	Standard Platform	3.75 / 4.2 / 5	MED257

# Lista de compatibilidad

Implant brand	System/model	Connection	Platform Type	Implant Diameter Ø	MED Nr.
ARDS	Premium	Internal Hex Connection	Standard Platform	3.75 / 4.2 / 5	MED257
ARDS	Classic		3.75	3.3/3.75/4.2/5/6	MED257
Astra	OsseoSpeed	Conical Connection	Wide Platform	4.5 / 5	MED332
Astra	Lilac		Lilac Wide	4.5/5	MED332
Avinent	Ocean	Internal Hex Connection		3.5 / 4	MED257
Avinent	Coral	Internal Hex Connection		3.5 / 4	MED257
Bioconcept	BV Tapered Bone Level		Narrow Platform	3.5	MED318
Bioconcept	BV Tapered Bone Level		Regular Platform	4/4.5/5	MED325
Bioconcept	BC Tissue Level Standard		Regular Platform	3.3/4.1/4.8	MED363
Bioconcept	BC Tissue Level Standard		Regular Platform	4.8	MED363
Bioconcept	BC Tissue Level Tapered		Regular Platform	4.8	MED363
Biogenesis	Aticon		Blue	3.5/4/4.5/5	MED332
Biogenesis	Iticon		4.8	3.5/4.1/4.8	MED363
BioHorizons	Tapered Internal plus	Internal Hex Connection		3.5	MED257
BioHorizons	Tapered Internal plus	Internal Hex Connection		4.5	MED387
Biohorizons	Tapered Internal		4.5 (Green)	4.6	MED349
Biohorizons	Tapered Plus		4.5 (Green)	5.8	MED349
Biohorizons	Mount-free Tapered		4.5 (Green)	4.6	MED349
Biohorizons	Tapered Pro		4.5 (Green)	5.2	MED349
Biohorizons	Tapered Short		4.5 (Green)	5.8	MED349
Biohorizons	Tapered Tissue Level		4.5 (Green)	4.6	MED349
Bioner	Hikelt	Internal Hex Connection	Regular Platform	3.8	MED257
Bioner	Hikelt	Internal Hex Connection	Wide Platform	4.7	MED387
Bioner	Hikelt		4.9	4.7	MED349
Biotec	SPR	Internal Hex Connection	Standard Platform	3.3 / 3.75 / 4.2 / 5 / 6	MED257
Biotec	SPTT	Internal Hex Connection	Standard Platform	4.2 / 5 / 6	MED257
Biotec	CIM	Internal Hex Connection	Standard Platform	3.3 / 3.75 / 4.2 / 5 / 6	MED257

# Lista de compatibilidad

Implant brand	System/model	Connection	Platform Type	Implant Diameter Ø	MED Nr.
Blue Sky Bio	Tri-Lobe (like Nobel)	Tri lobe	Narrow Platform	3.5	MED271
Blue Sky Bio	Tri-Lobe (like Nobel)	Tri lobe	Regular Platform	4.3	MED295
Blue Sky Bio	Tri-Lobe (like Nobel)	Tri lobe	Wide Platform	5	MED301
Blue Sky Bio	Conus	Conical Connection	Wide Platform	4.5 / 5	MED332
BTK	Klassic / Konic		3.5 IR	3.25/4	MED257
BTK	IS +		DR	3.3/3.7/4.1/4.8/6	MED318
BTK	BT Safe BL		KR	3.3/3.7/4.1/4.8	MED318
Camlog	Screw-Line	Camlog Connection-3 slots		5	MED639
Camlog	Screw Line	Camlog Connection-3 slots		3.8	MED622
Camlog	Screw Line	Camlog Connection-3 slots		4.3	MED615
Cortex	Dynamic	Internal Hex Connection	Standard Platform	3.8 / 4.2 / 5	MED257
Cortex	Classix	Internal Hex Connection	Standard Platform	3.8 / 4.2 / 5	MED257
Cortex	Saturn	Internal Hex Connection	Standard Platform	3.8 / 4.2 / 5	MED257
Cortex	Conical Dynamix		Wide Platform	5.0/6.0	MED332
Cortex	Conical Classix		Wide Platform	5.0/6.0	MED332
C-Tech	Esthetic Line	Conical Connection	Standard Platform	3.3	MED660
DENTAL TECH	Implogic		4.5 (Blue)	4.5	MED349
DENTEGRIS	SLS-Straight		4.5	4.5	MED349
DENTEGRIS	Sinus-Lift		4.5	4.5	MED349
DENTEGRIS	S&T Implants		4.5	4.5	MED349
Dentis	s-Clean Save		5.5/6	5.5/6	MED325
Dentis	SQ-SL		Regular Platform	4/4.5/5	MED325
Dentis	SQ-SL		Wide Platform	6/7/8	MED325
Dentis	i-Clean Tapered		Narrow Mini	3.7	MED363
Dentis	i-Clean Tapered		Regular Platform	4.1/4.3	MED363
Dentis	i-Clean Tapered		Wide Platform	4.8	MED363
Dentis	i-Clean Straight		4.8	4.1/4.8	MED363

# Lista de compatibilidad

Implant brand	System/model	Connection	Platform Type	Implant Diameter Ø	MED Nr.
Dentium	Superline	Internal Hex Connection	All platforms	3.6 / 4.0 / 4.5 / 5 / 6 / 7	MED325
Dentium	SimpleLine II		4.8	3.8/4.3	MED363
Dentsply	Ankylos C/X	Tapered	Regular Platform	3.5	MED271
DIO	UF II		Regular Platform	3.8/4.4/5/5/5.5	MED325
DMI DENTAL SUPPLY	DCI/DSI		3.75	3.3/3.5/3.75/4.2/5/6	MED257
Easy Implants	Master C		4.5 (Lilas)	4.5	MED325
Easy Implants	Master C		5 (Lilas)	5	MED325
Easy Implants	Master S		3.75 (Lilas)	3.75	MED325
Easy Implants	Master S		4.25 (Lilas)	4.25	MED325
Easy Implants	Master S		4.75 (Lilas)	4.75	MED325
Easy Implants	Master L		3.3 (Lilas)	3.3	MED325
Easy Implants	Master L		3.75 (Lilas)	3.75	MED325
Easy Implants	Master L		4.25 (Lilas)	4.25	MED325
Easy Implants	Master L		4.75 (Lilas)	4.75	MED325
Eckermann	Winner		3.5	3/3.5/4	MED257
Eckermann	Winner		4.5	4/4.5/5	MED349
Euroteknika	Aesthetica		4.8	4.1	MED363
Euroteknika	Ibone G		Regular Platform	4.8/5.5	MED363
F&B Implant (Fit & Brilliant)	FA Submerged Fixture		Narrow Platform	3.9	MED325
F&B Implant (Fit & Brilliant)	FA Submerged Fixture		Regular Platform	4.1/4.4	MED325
F&B Implant (Fit & Brilliant)	FA Submerged Fixture		Wide Platform	4.8	MED325
F&B Implant (Fit & Brilliant)	FA Submerged Fixture		Ultra-Wide Platform	5.3/5.8/6.3/6.8	MED325
Glidewell	Inclusive	Internal Hex Connection	Standard Platform	3.3 / 3.7 / 4.1	MED257
Glidewell	Inclusive	Internal Hex Connection	Regular Platform	4.5	MED387
GMI (ilerimplant)	Frontier		Wide Platform 4.3	4.75/5.75	MED349
GMI Ilerimplant	Frontier		Regular Platform 3.3	3.3/3.75/4.25	MED257
GT Medical	Best Fit Hexágono Interno		Wide Platform	3.7/4.1/4.3/4.8	MED332

# Lista de compatibilidad

Implant brand	System/model	Connection	Platform Type	Implant Diameter Ø	MED Nr.
GT Medical	Best Fit Octógono Interno		Regular Platform	3.7/4.3/4.8	MED363
Hahn Implant - Glidewell	Hahn Tapered Implant		5	5	MED264
Hahn Implant (Glidewell)	Hahn Tapered Implant		5	5	MED356
Hiossen	ET III SA/BA, ET II SA	Conical Connection	Regular Platform H(2.1)	3.5	MED318
Hiossen	ET III SA/BA, ET II SA	Conical Connection	Wide Platform H(2.5)	4 / 4.5 / 5	MED325
HI-TEC	Tapered Self Thread		3.5	3.3/3.75	MED257
HI-TEC	Logic Plus		3.7	3.5	MED257
HI-TEC	Logic Plus		3.9	4.3	MED257
HI-TEC	Tapered Self Thread		4.5	4.2/5	MED349
IBS	Magic FC		3.8	4/4.5/5/5.5/6/6.5	MED325
IBS	N.R. Fix		3.8	3/3.5	MED325
IDO Implants	I Do Implant		Universal	3.8/4/4.5/5/5.5/6/7	MED325
Implanet	Neo	Internal Hex Connection	Standard Platform	3.5 / 3.75 / 4.2 / 5 / 6	MED257
Implant Direct	InterActive	Conical Connection		4.3 / 5	MED264
Implant Direct	Tri-Lobe (like Nobel)	Tri Lobe	Regular Platform	4	MED288
Implant Direct	Replant / Replus / ReActive	Tri Lobe	Narrow Platform	3.5 / 3.7 / 4.2	MED271
Implant Direct	Replant / Replus / ReActive	Tri Lobe	Regular Platform	4.3	MED295
Implant Direct	Replant / Replus / ReActive	Tri Lobe	Wide Platform	5	MED301
Implant Direct	Legacy 1 / 2 / 3 / 4	Internal Hex Connection	Standard Platform	3.7 / 4.2	MED257
Implant Direct	Legacy 1 / 2 / 3 / 4	Internal Hex Connection	Standard Platform	4.7 / 5.2	MED387
Implant Direct	SwishActive		3.4	4.1/4.8	MED264
Implant Direct	Simply Iconic		Gold 3.4	4.7/5.2/5.7	MED264
Implant Direct	Legacy		4.5	4.7/5.2	MED349
Implant Direct	SwishActive		3.4	4.1/4.8	MED356
Implant Direct	Interactive		3.4	4.3/5	MED356
Implant Direct	Simply Iconic		Gold 3.4	4.7/5.2/5.7	MED356
Implant Direct	Swishplant / Swishplus		4.8	4.1/4.8	MED363

# Lista de compatibilidad

Implant brand	System/model	Connection	Platform Type	Implant Diameter Ø	MED Nr.
Implant Swiss	Bone Level		3.7	3.7	MED325
Implant Swiss	Bone Level		4.3	4.3	MED325
Implant Swiss	Bone Level		4.8	4.8	MED325
Implant Swiss	Bone Level		5.5	5.5	MED325
Intra-Lock	IntraHex		3.5	3.75/4	MED257
Intra-Lock	IntraHex		4.5	4.75	MED349
J Dental Care	JD Evolution	Internal Hex Connection	Standard Platform	3.7 / 4.3 / 5 / 6	MED257
J Dental Care	Plus	Internal Hex Connection	Standard Platform	3.7 / 4.3 / 5 / 6	MED257
J Dental Care	JD Pterygo		4	4	MED257
J Dental Care	JD ICON		3.9 - 4 - 4.7	3.9 / 4.3 / 5	MED264
J Dental Care	JD ICON		3.9 - 4 - 4.7	3.9 / 4.3 / 5	MED356
J Dental Care	JD ICON Plus		3.7 - 4 - 4.8	3.7 / 4.3 / 5	MED575
J Dental Care	JD ICON Plus T		3.5	3.5 / 4 / 4.5	MED575
Kyocera	Finesia	Internal Hex Connection	Narrow Platform	3.5	MED271
LEADER	Tixos Internal Hex		3.5	3.3	MED257
LEADER	Tixos Internal Hex		4	3.75	MED257
Megagen	AnyRidge	Internal Hex Connection	Regular Platform	3.5 / 4 / 4.5 / 5 / 5.5 / 6	MED677
Megagen	AnyOne Internal		Regular Platform	3.5/4/4.5/5/6/7	MED325
Megagen	ExFeel		Small	3.5	MED363
Megagen	ExFeel		Regular Platform	4.1	MED363
Megagen	ExFeel		Wide Platform	4.8/5	MED363
Meisinger	Tissue Level	Internal Hex Connection	Regular Platform	3.3. 3.75. 4.1. 4.8	MED363
Meisinger	Tissue Level	Internal Hex Connection	Wide Platform	4.8	MED363
Microdent	Ektos		3.5	3.7/4.2	MED257
MIS	Seven	Internal Hex Connection	Standard Platform	3.75 / 4.2	MED257
MIS	Seven	Internal Hex Connection	Wide Platform	5 / 6	MED349
MIS	M4	Internal Hex Connection	Standard Platform	3.75 / 4.2 / 5	MED257

# Lista de compatibilidad

Implant brand	System/model	Connection	Platform Type	Implant Diameter Ø	MED Nr.
MIS	C1	Conical Connection	Standard Platform	3.75 / 4.2	MED394
MIS	C1	Conical Connection	Wide Platform	5	MED400
MIS	M4		Wide Platform	5.0/6.0	MED349
MIS	V3		Standard Platform	3.9/4.3/5	MED394
MPI	Privilege CM		Wide Platform	5	MED332
MPI	Excellence CM		Wide Platform	5	MED332
MSDI	Alef System	Internal Hex Connection	Standard Platform	3.5. 3.75. 4.2. 5. 6	MED257
Neobiotech	IS Implant System		S-Narrow Platform	3.2	MED318
NeoBiotech	IS Implant System		Narrow Platform	3.5	MED325
NeoBiotech	IS Implant System		Regular Platform	4	MED325
NeoBiotech	IS Implant System		Regular Platform	4.5	MED325
NeoBiotech	IS Implant System		Wide Platform	5	MED325
<a href="#">Nobel</a>	NobelActive TiUltra	Internal Hex	Regular Platform	4.3 / 5	MED264
<a href="#">Nobel</a>	NobelParallel CC TiUltra	Conical Connection	Regular Platform	4.3 / 5	MED356
Nobel	NobelActive	Internal Hex	Regular Platform	4.3 / 5	MED264
Nobel	NobelParallel CC	Conical Connection	Regular Platform	5	MED356
Nobel	NobelReplace CC	Conical Connection	Narrow Platform	3.5	MED271
Nobel	NobelReplace CC	Conical Connection	Regular Platform	4	MED288
Nobel	NobelReplace CC	Conical Connection	Regular Platform	4.3	MED295
Nobel	NobelReplace CC	Conical Connection	Regular Platform	5	MED301
Nobel	NobelParallel CC	Conical Connection	Regular Platform	4.3	MED356
Noris Medical	Onyx	Internal Hex Connection	Standard Platform	3.3 / 3.75 / 4.2 / 5 / 6	MED257
Noris Medical	Tuff Pro	Internal Hex Connection	Standard Platform	3.3 / 3.75 / 4.2 / 5 / 6	MED257
Noris Medical	Tuff TT	Internal Hex Connection	Standard Platform	4.2 / 5 / 6	MED257
Noris Medical	Tuff	Internal Hex Connection	Standard Platform	3.3 / 3.75 / 4.2 / 5 / 6	MED257
Normon	Normoimplant HI		3.5	3.75/4.25/4.75	MED257
NOVA IMPLANTS	PSI/PCI		3.75	3.3/3.75/4.2/5/6	MED257

# Lista de compatibilidad

Implant brand	System/model	Connection	Platform Type	Implant Diameter Ø	MED Nr.
Osstem	TS III SA	Internal Hex Connection	Regular Platform	4 / 4.5 / 5	MED325
Osstem	TS III SA	Internal Hex Connection	Ultra Wide Platform	6 / 7	MED325
OSSTEM IMPLANT	TS		Narrow Mini 3.5	3.5	MED318
Osteoplus	Shi		3.5	3.3 / 3.75 / 4.2	MED257
OXY	K1 Line		Regular Platform	3.5/4/4.5/5	MED677
OXY	K1 Line		Wide Platform	5.5/6/6.5	MED677
OXY	PSK Line		Regular Platform	3.5/4/4.5/5	MED677
OXY	MD Line KONE		Regular Platform	3.75/4.25/5	MED677
Paltop	Advanced Classic	Internal Hex Connection	Standard Platform	3.75 / 4.2 / 5 / 6	MED257
Paltop	Advanced +	Internal Hex Connection	Standard Platform	3.75 / 4.2 / 5 / 6	MED257
Paltop	Dynamic	Internal Hex Connection	Standard Platform	3.75 / 4.2 / 5 / 6	MED257
Paltop	DIVA	Internal Hex Connection	Standard Platform	3.75 / 4.2 / 5	MED257
Paltop	Conical Active		Standard Platform	3.25/3.75/4.2/5	MED318
Paltop	Advanced classic		Wide Platform (purple)	6	MED349
Paltop	Advanced +		Wide Platform (purple)	6	MED349
Paltop	Dynamic		Wide Platform (purple)	6	MED349
Paltop	PAI		Wide Platform (Purple)	6	MED349
Proclinic	Cilindrico Interno/Cónico		3.5	3.3/3.75/4.25/5/3.5/4/5	MED257
Proclinic	SP Octa		4.8	3.3/4.1/4.8	MED363
Radhex	PHI		3.5	3.75	MED257
Radhex	PHI		4.5	4.5/5	MED349
Reflect	Tapered Screw		3.5	3.5	MED257
Reflect	Rapid		Regular Platform	4.3/5/5.5	MED264
Reflect	Recover		Regular Platform	4.3/5	MED264
Reflect	Aspire		Lilac (Ancha)	5	MED332
Reflect	Tapered Screw		4.5	4.1/4.7	MED349
Reflect	Rapid		Regular Platform	4.3/5/5.5	MED356

# Lista de compatibilidad

Implant brand	System/model	Connection	Platform Type	Implant Diameter Ø	MED Nr.
Reflect	Recover		Regular Platform	4.3/5	MED356
SEWON MEDIX	IH2 SLA SYSTEM		Narrow Mini	3.5	MED318
SEWON MEDIX	IH2 RBM SYSTEM		Narrow Mini	3.5	MED325
Sewon Medix	IH2 SLA SYSTEM		Regular Platform	3.5/4/4.5/5	MED325
Sewon Medix	IH2 RBM SYSTEM		Regular Platform	3.5/4/4.5/5	MED325
SGS Dental	P1	Internal Hex Connection	Standard Platform	3.75 / 4.2 / 5 / 6	MED257
SGS Dental	P7	Internal Hex Connection	Standard Platform	3.75 / 4.2 / 5 / 6	MED257
Signo Vinces	Duocon		CM4.6	4.6/5.5	MED332
Southern Implants	IT - Cylindrical &Tapered	IT4.8, IT6.5		3.3. 4.1. 4.9	MED363
Southern Implants	Internal Hex - Tapered	IM-M		3.75. 4.20. 5.00	MED257
Southern Implants	Provata - Tapered	PRO4, PRO5		4.0 / 5.0 / 6.0	MED257
Southern Implants	Provata - Tapered	PRO3		3.3	MED257
Southern Implants	Internal Hex		Universal	3.75/4.2/5	MED257
Southern Implants	Deep Conical		4.1	5	MED332
Southern Implants	IT Connection		4.8	3.3/4/4.1/4.9/5	MED363
Sterngold	PUR	Internal Hex Connection	Narrow Platform	3.5 / 4.3	MED257
Sterngold	PUR	Internal Hex Connection	Standard Platform	5.0 / 6.0	MED387
Sterngold	TRU	Conical Connection	Narrow Platform	3.5 / 4.3 / 5 / 6	MED264
Straumann	Tissue Level	Standard Implant	Standard Implant	3.3 / 4.1 / 4.8	MED363
TICARE	MG Inhex		3.8 Maxi	5	MED332
Titanium Fix	b-fix		Larga	4.5/5	MED332
Tov Implants	Maer	Internal Hex Connection	Standard Platform	3.5 / 3.75 / 4.2 / 5 / 6	MED257
Tov Implants	Ragil	Internal Hex Connection	Standard Platform	3.5 / 3.75 / 4.2 / 5 / 6	MED257
TREE-OSS	Simple		3.75 Amarillo (Yellow)	3.3/3.75/5	MED257
Tree-oss	Anatomic/HS		3.5 Rose	3.5	MED271
Tree-oss	Anatomic/HS		4.3 Yellow	4.3	MED295
Tree-oss	Anatomic/HS		5 Blue	5	MED301

# Lista de compatibilidad

Implant brand	System/model	Connection	Platform Type	Implant Diameter Ø	MED Nr.
TRI Dental Impants	TRI-Vent (Performance)	Internal Hex Connection	Standard Platform	3.75 / 4.1 / 4.7	MED257
TRI DENTAL implants	TRI-Vent		3.5	3.75/4.1/4.7	MED257
UFIT	Gt2		Regular Platform	4/4.5	MED332
UFIT	Gt2		Wide Platform	5	MED332
UFIT	Gt2		Ultra-Wide Platform	5.5/6/6.5/7	MED332
UFIT	Nt2		Regular Platform	4/4.5	MED332
UFIT	Nt2		Wide Platform	5	MED332
UFIT	Nt2		Ultra-Wide Platform	5.5/6/6.5/7	MED332
VULKAN IMPLANTS	IN-Hex		3.75	3.3/3.75/4.2/5	MED257
XGATE	Pure & Porous	Internal Hex Connection	Standard Platform	3.3 / 3.75 / 4.2 / 5	MED257
Yes Implant	S-SYSTEM		Narrow Platform	3.3/3.5	MED325
Yes Implant	S-SYSTEM		Regular Platform	4/4.5	MED325
Yes Implant	S-SYSTEM		Wide Platform	5/5.5	MED325
Ziacom	ZINIC		Regular Platform 3.5	3.7/4/4.3	MED257
Ziacom	ZINIC Shorty		Regular Platform 3.5	4.75	MED257
Zimmer	TSV	Internal Hex Connection	Standard Platform	3.5 / 4.1 / 4.3 / 5.1	MED257
Zimmer	MT	Conical Connection		4.1 / 5.1	MED417
Zimmer	Screw-Vent		3.5	3.7/4.1	MED257
Zimmer	TSX Implant		3.5	3.7/4.1/4.7	MED257
Zimmer	Screw-Vent		4.5	4.7	MED349
Zimmer	TSX Implant		4.5	5.4/6	MED349
Zimmer	Swiss-Plus		4.8	3.7/4.1/4.8	MED363



**Magdent Ltd.**

*Address*

B.S.R. 3, Kinneret St 5, Bnei Brak, Israel

*Phone*

+972-524-329-329

*Email*

[service@magdentmed.com](mailto:service@magdentmed.com)

*Website*

[www.magdentmed.com](http://www.magdentmed.com)